



WATREX PRICE LIST 2010

AMINO ACID - AQC kit

P/N	Description	Price CZK
D2K4101926	AQC kit for aminoacids determination Kit includes: D2K4101927 - 5x1 ml derivatizing agents for aminoacids determination for 5x50 analysis solution (6-aminochinolyl-N-hydroxysukcinimidyl carbamate in acetonitrile) D2K4101928 Column WATREX 250x4 mm AminoAcid-AQC 5 µm D2K4101929 3 pcs of guard columns WATREX 20x4 mm AminoAcid-AQC 5 µm WI029012 WATREX integral column holder Manual for aminoacids derivatisation	22.545,00
D2K4101927	Derivatizing agents AQC 5x1 ml solution of 6-aminochinolyl-N-hydroxysukcinimidyl carbamate in acetonitrile for aminoacids derivatisation	7.150,00
D2K4101928	WATREX 250x4 mm C18-AQ 5 µm Column -silica sorbent with octadecyl groups, suitable for separation of AQC-aminoacids derivatives	11.725,00
D2K4101929	Guard column insert 20x4 mm C18-AQ 5 µm set of 3 pcs	2.150,00
WI029012	Integral guard column holder 20 mm	2.700,00